



Testing Adapters for SPECTROMAXx

Different Adapters for Different Requirements



Small Parts Adapter Kit
(Part Nr. 75260151)


The SPECTROMAXx covers a wide range of analysis requirements. One major challenge relates to the various shapes and sizes of the samples and the optimized position on the spark stand. The new adapter kit is focused on ease of use and flexibility.




Spark Stand Plate with 8 mm electrode hole



Adapter for Wires with a diameter of 1-3 mm



Adapter for Foils or small pieces



Adapter for Wires/Screws or small parts with a diameter of 3-10 mm
(Part Nr. 75060178)

www.spectro.com

AMETEK[®]
MATERIALS ANALYSIS DIVISION

GERMANY

SPECTRO Analytical Instruments GmbH
Boschstrasse 10
D-47533 Kleve
Tel: +49.2821.8922102
Fax: +49.2821.8922202
spectro.sales@ametek.com

U.S.A.

SPECTRO Analytical Instruments Inc.
91 McKee Drive
Mahwah, NJ 07430
Tel: +1.800.548.5809
+1.201.642.3000
Fax: +1.201.642.3091
spectro-usa.sales@ametek.com

Hong Kong (Asia-Pacific)

SPECTRO Analytical Instruments
(Asia-Pacific) Ltd.
Unit 1603, 16/F, Tower III Enterprise Sq.
No. 9 Sheung Yuet Road
Kowloon Bay, Kowloon
Tel: +852.2976.9162
Fax: +852.2976.9132
spectro-ap.sales@ametek.com

Subsidiaries: **CHINA:** Tel. +86.10.8526.2111, Fax +86.10.8526.2141, spectro-china.sales@ametek.com, **FRANCE:** Tel. +33.1.3068.8970, Fax +33.1.3068.8999, spectro-france.sales@ametek.com, **GREAT BRITAIN:** Tel. +44.1162.462.950, Fax +44.1162.740.160, spectro-uk.sales@ametek.com, **INDIA:** Tel. +91.22.6196.8200, Fax +91.22.2836.3613, sales.spectroindia@ametek.com, **ITALY:** Tel. +39.02.94693.1, Fax +39.02.94693.650, spectro-italy.sales@ametek.com, **JAPAN:** Tel. +81.3.6809.2405, Fax +81.3.6809.2410, spectro-japan.sales@ametek.co.jp, **SOUTH AFRICA:** Tel. +27.11.979.4241, Fax +27.11.979.3564, spectro-za.sales@ametek.com, **SWEDEN:** Tel. +46.8.5190.6031, Fax +46.8.5190.6034, spectro-nordic.sales@ametek.com.

SPECTRO operates worldwide and is present in more than 50 countries. For SPECTRO near you please visit www.spectro.com/worldwide.

© 2012 SPECTRO. Subject to modifications • H-12, Rev.0. Photos: SPECTRO, Corbis, Getty Images, iStockphoto.

Registered trademarks of SPECTRO Analytical Instruments GmbH **SPECTRO**: USA (3,645,267); EU (005673694); "SPECTRO": EU (009693763); iCAL: USA (3,189,726), EU (003131919);